


BEAM APPARATUS

VDAS® SM1004

Examines the deflection and forces on different types of beams for a wide range of supports and loads; also demonstrates Young's modulus.



SCREENSHOT OF THE OPTIONAL VDAS® SOFTWARE

- Includes textbook with full theory
- Simply supported and cantilever beam tests with up to four supports with any loading
- Three load cells with digital indicators measure reaction forces or act as rigid sinking supports
- Precision digital indicators for accurate deflection measurements

LEARNING OUTCOMES:

- Verification of the bending equation
- Determination of flexural rigidity and elastic modulus (Young's modulus)
- Verification of static equilibrium
- Deflection of beams on two simple supports with point loads
- Reciprocal properties for loads and deflection
- Simple and propped cantilevers with any loading
- Continuous beams – statically indeterminate cases for simply supported beams and cantilevers on more than two supports with any loading (including measurement of unknown reactions)
- Simply supported and cantilever beams with sinking supports

With the SM1004a Specimen Beams, these additional experiments can be done:

- The effects of material and section shape on flexural rigidity
- Bending characteristics of a brass/steel compound beam, with and without shearing connection between the two layers
- Equivalent sections – characteristics of a metal-faced wooden beam
- Deflections on a non-uniform (tapered) beam or cantilever

CONTINUED ON NEXT PAGE

BEAM APPARATUS (SM1004) CONTINUED FROM PREVIOUS PAGE

The apparatus consists of an upper cross-member carrying graduated scales, and two lower members bolted to T-legs to form a rigid assembly. The three load cells and cantilever-support pillar slide along the lower members and can be clamped firmly in any position. The load cells have direct digital readout and each is fitted with a hardened steel knife edge which can be adjusted to set the initial level, or to simulate a sinking support. Locking pins can convert each load cell to a rigid support when required. The cantilever support is a rigid pillar with a sturdy clamping arrangement to hold the beams when built-in end conditions are required. Four weight hangers and a set of weights are supplied to apply static loads.

RECOMMENDED ANCILLARIES:


- Versatile Data Acquisition System – Bench-mounted version (VDAS-B) 299
- Additional Specimen Beams (SM1004a)

ALTERNATIVE PRODUCTS:


- Deflection of Beams and Cantilevers Kit (ES4) 11
- Stiffness – Bending and Torsion (TE16) 152
- Beam and Leaf Spring (SM1000g) 165
- Deflection of Beams and Cantilevers (STR4) 202
- Continuous and Indeterminate Beams (STR13) 204

PLAY NOW - TECQUIPMENT YOUTUBE CHANNEL

Visit and subscribe to the TecQuipment YouTube channel for all the latest products showcase, installation and user videos. [YOUTUBE.COM/C/TECQUIPMENT](https://www.youtube.com/c/tecquipment)



GEAR TRAINS (ES13)




AIR CONDITIONING TRAINER (EG1501)




ROTATING FATIGUE MACHINE (SM1090)



HYDROLOGY AND RAINFALL (H313)



TRACE DAMPING MODELS (VDAS®)



SMALL ENGINE TEST SET (TD200)